

# If an outlet looks like this, chances are the Integrated Flowmeter is completely compatible



**Take a picture of the outlet and back body and send it to your Amvex Sales Manager to be sure.**

TRUSTED BRANDS OF OHIO MEDICAL

aeros<sup>®</sup>  
Division of Ohio Medical


Amvex Corporation  
Division of Ohio Medical

Ohio  
Medical Corporation™

Squire-Cogswell  
Division of Ohio Medical

# Not Compatible As Of Now

If you see any of these outlets they may not be compatible without a retrofit kit



Puritan Bennett



Allied/NCG/Chemtron



Medstar / Oxequip



Ohmeda  
Diamond 1

TRUSTED BRANDS OF OHIO MEDICAL

# Information We Need:

- **Facility Name**
- **End User Names – Director of RT, Maintenance/Plant Ops, Other Contacts**
- **Wall Outlet Manufacturer Type (eg: Beacon)**
- **Outlet Connection Type (eg: Ohmeda)**
- **New Construction or Retrofit opportunity?**
- **Volume + Hospital Department conducting evaluation**
- **Take a picture of the outlet face plate then remove the plate and take a picture of the back assembly**

TRUSTED BRANDS OF OHIO MEDICAL

**aeros**<sup>®</sup>  
Division of Ohio Medical

**A** Corporation  
**AMVEX**  
Division of Ohio Medical

**Ohio**  
Medical Corporation™

**Squire-Cogswell**  
Division of Ohio Medical